



Application

DC Tests (I-V/C-V, PIV) , AC tests, etc.

Features

1. Small footprint (6cm*4cm*5cm).
2. Very cost-effective micropositioner.
3. Applied for pads down to 50 μ m*50 μ m.
4. Standard products can be shipped within 2 weeks.

Specifications

Accuracy	10 μ m or 635 μ m/rev
X-Y-Z Travel Range	12.5mm*12.5mm*12.5mm
TPI	40
Mounting	Magnetic (Without On/Off Switch)
Dimensions	60mm*40mm*50mm (W*D*H)

*Please contact for detailed information of this model.

*Please note that specifications and design are subject to change without notice.