



### Application

DC Tests (I-V/C-V, PIV) , AC tests, etc.

### Features

1. Small footprint (8cm\*6cm\*5cm).
2. Very cost-effective micropositioner.
3. Applied for pads down to 30µm\*30µm.

### Specifications

Accuracy	5µm or 500µm/rev
X-Y-Z Travel Range	12.5mm*12.5mm*12.5mm
TPI	50
Mounting	Magnetic-ON/OFF
Dimensions	80mm*60mm*50mm (W*D*H)

\*Please contact for detailed information of this model.

\*Please note that specifications and design are subject to change without notice.